

# Total Chlorine Testing

## Pitfalls And Remedies

*Presented By*



NATIONAL  
DIALYSIS  
ACCREDITATION  
COMMISSION<sup>®</sup>

# What Is It & Why Is It In Our Water?

Chlorine (e.g., bleach) is a disinfectant; chloramine is a form of chlorine that stays active longer

Cities add chlorine and chloramines to their water distribution systems to protect our drinking water from bacterial growth

# Why An Accurate Test Is Important: Patient Consequences

The blood of hemodialysis patients  
is exposed to a LOT of water

24 – 50 gallons each treatment

Patient injury or death can happen  
when water used to make dialysate  
contains even a low level of total  
chlorine

# Why An Accurate Test Is Important: Dialysis Facility Consequences



- Condition-level deficiencies result if:
  - You do not follow the test directions EXACTLY
  - You do not properly document the test for total chlorine
  - You use expired testing strips
  - You use a strip or reagent that is not sensitive to the required levels
- At an initial survey, Condition-level deficiencies = **failed survey**
- At any survey, Condition-level deficiencies = a repeat visit

# Total Chlorine Testing

- By using the cut-off of the maximum allowable level for chloramine (<0.1 ppm)
  - ☐ One test checks for both chlorine AND chloramine
  - ☐ **You protect patients from both**
- Total chlorine test methods:
  - ☐ Dip-and-read test strips
  - ☐ “DPD” based kits (tablets or powders-rarely used now)
  - ☐ On-line monitoring
- Accurate results depend on following the directions EXACTLY



Location

Test Sensitivity

Sample

Timing

Documentation

## Pitfalls Of Total Chlorine Testing

# Pitfall: Location

Where should the sample be taken?

- Routine samples should be taken from a port AFTER the first carbon filter


What if that sample is  $>0.1$ ?

- If that sample shows a level  $>0.1$  mg/L, sample from a port after the second carbon filter

Remedy: Label collection ports and follow facility policy

# Pitfall: Test Sensitivity

- Total chlorine maximum “safe” level is less than or equal to 0.1 mg/L
- Test used must measure levels LOWER than 0.1 mg/L
- Zero  $\neq$  zero:
  - Tests only measure to “less than” the lowest number on the color comparison chart
    - Example: if color chart shows 0.5 and then 0, the test only measures to “less than 0.5”
  - No test measures to “zero”
- Warning: tests for residual chlorine usually measure only to 0.5mg/L



Remedy: Do not store residual test strips in total chlorine testing area

Use only methods that test to less than 0.1 mg/L




# Pitfall: Test Sensitivity

Expiration date:


- Testing materials (test strips, powders, tablets) out-of-date

Remedy: Always check packaging for expiration date

- Discard expired strips
- Alert management to reorder before need is critical



# Pitfall: Sample



# Pitfall: Timing

Testing frequency:

- At the start of each treatment day prior to first patient treatment
- Prior to each shift or every 4 hours if there are no set patient shifts

Remedy: follow facility policy for test frequency

# Pitfall: Timing



Test water that has spent the normal time in the carbon filter



Removal of total chlorine by carbon requires 10 minutes



When the RO is not running, water remains in carbon tank longer

Testing a water sample that was in the carbon tank longer = falsely lower value, potential for patient harm



Remedy: Ensure RO has been running at least 15 minutes prior to collecting sample

# Pitfall: Timing

When performing the test:

- USE A TIMING DEVICE (watch with second hand or clock **or timer**)
- Each step has a specific time requirement
  - ☐ Test is only accurate if **time** requirements are met
- Meet each time requirement EXACTLY
  - ☐ Examples:
    - ☐ Expose strip to sample for **60** seconds
    - ☐ Wait **20** seconds before reading

Remedy: Invest in a clock with a second hand **or timer** for the water room

Educate all who do testing to the importance of exact timing

# Online Total Chlorine Monitoring



- Monitoring is continuous, so no particular testing times are required
- Manufacturer requires a manual test of total chlorine before the first patient treatment each day to verify the accuracy of the on-line system
- Facility policy may require documentation of the monitor values



# Timing

Comparing the test to the color chart

- Must use a **white** background
- Must be done immediately after test completed
  - Color can change when reading is delayed
- If second person verification is required, second person should be present when the test is performed
  - Take sample and testing equipment to the treatment floor if necessary
- Test must always be compared to the color chart

Remedy: Keep color chart handy

Educate all who do testing to use color chart **promptly**

# Pitfall: Documentation



- Don't document the test until it is completed
- Be sure to document accurately (**correct:**  $<0.1$ , **not correct:**  $>0.1$  or  $1.0$ )
- Careful in documenting time and date

Remedy: Educate all who do testing to the need for careful documentation

Audit: Both practice and documentation



# Avoiding Pitfalls

*“There is a right way, staff know the right way, and staff perform the task the right way.”*

- Facility policy should:
  - ✓ Follow the manufacturer’s guidance
  - ✓ Provide clear step-by-step direction
- Staff should be trained on the “whys” of each step in total chlorine testing
- Competency testing is required for each staff member assigned the task
- Audits of practice provide opportunity to identify and correct any issues
- It’s okay to use the posted procedure for reference: may help with nervousness

# Supervision

- The registered nurse in charge needs to know the basics of water treatment and ensure that assigned staff members are accurately performing the tests and recording the results
- The Medical Director and Nurse Manager must ensure all staff assigned water treatment duties are qualified by training and have current documentation of competency



# Summary

- Safe water for treatment is critical to quality care and safety for every dialysis patient
- Staff members who understand the importance of accurate total chlorine testing will keep patients safe and clinics survey ready



# References

- Centers for Medicare & Medicaid Services (CMS). (2008a). *Conditions for Coverage for end-stage renal disease facilities; Final Rule*. <http://www.cms.gov/regulations-and-guidance/legislation/cfcsandcops/downloads/esrdfinalrule0415.Pdf>
- Centers for Medicare and Medicaid Services (CMS). (2008b). *ESRD Surveyor Training; Interpretive Guidance*. <https://www.cms.gov/medicare/provider-enrollment-and-certification/guidanceforlawsandregulations/downloads/esrdpgmguidance.Pdf>
- Hellebrand, A., Allen, D., & Hoffman, M. (2017). Hemodialysis. In S.M. Bodin (ed.), *Contemporary Nephrology Nursing* (3<sup>rd</sup> ed., Pp. 153-207). American Nephrology Nurses Association.
- International standards organization (ISO). (2019). *ISO 23500-3: Preparation and quality management of fluids for haemodialysis and related therapies- Part3: Water for haemodialysis and related therapies*. <https://www.iso.org/standard/67612.Html>
- Medical Education Institute (MEI). (2018). *Core curriculum for the dialysis technician* (6<sup>th</sup> ed.). The Medical Education Institute, Inc.
- National Dialysis Accreditation Commission (NDAC). (2018) *Standards and Interpretive Guidance*. Available on request from [info@ndacommission.com](mailto:info@ndacommission.com).

# References

- Payne, G.M. (2020). Hemodialysis: water treatment. In C.S. Counts (ed.), *Core Curriculum for Nephrology Nursing* (7<sup>th</sup> ed., Pp. 988-1001). American Nephrology Nurses Association.
- Payne, G.M., & Curtis, J. (2018). Water treatment for hemodialysis: what you must know to keep patients safe. *Nephrology Nursing Journal*, 45(2), 141-168.
- U.S. Food and Drug Administration (FDA). (1997). *Guidance for the content of premarket notifications for water purification components and systems for hemodialysis*. <https://www.fda.gov/media/72713/download>
- Payton, J. (2021). *Clinical Consult: What is the importance of testing to total chlorine in water used for hemodialysis?* *Nephrology Nursing Journal*, 48(2), 189-192.